

Datasheet: PFP003

BATCH NUMBER 164095

Description:	RECOMBINANT CAT SERUM AMYLOID A
Name:	SERUM AMYLOID A
Other names:	SAA
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			
Functional Assays			■	

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Cat
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 1.0 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified by affinity and reversed-phase chromatography
Preservative Stabilisers	None present
Approx. Protein Concentrations	0.1 mg/ml after reconstitution
External Database Links	UniProt:

Product Information **Recombinant Cat serum amyloid A (SAA)** is a 121 aa recombinant protein produced in *E.coli*, which includes a proprietary 10 aa tag at the N-terminus.

Serum amyloid A proteins are a family of acute phase proteins, produced primarily by the liver, highly conserved in vertebrate species, and used as markers of inflammation.

Feline SAA, as with serum amyloid A in other species, is present in different isoforms, this recombinant protein corresponds to isoform 1 (SAA1). Feline SAA is composed of 111 amino acids and differs from human SAA by having an 8 aa sequence insert in the central region of the sequence.

With similar biochemical properties to the native feline SAA, this reagent may be used as a standard in immunoassays for measurement of levels of cat serum amyloid A.

Protein Molecular Weight	13.8 kDa
---------------------------------	----------

Purity	>95%
---------------	------

Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
----------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	------------------------------------------------------------

Health And Safety Information	Material Safety Datasheet documentation #10268 available at: https://www.bio-rad-antibodies.com/SDS/PFP00310268
--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Regulatory	For research purposes only
-------------------	----------------------------

North & South America Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M412105:221109'

Printed on 18 Nov 2024