

Datasheet: PHP046

BATCH NUMBER 155864

Description:	RECOMBINANT HUMAN INTERLEUKIN-7
Name:	IL-7
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4 ng/well
Western Blotting	▪			1.5 - 3.0 ng/lane
Functional Assays	▪			0.1ng/ml - 10ng/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the protein for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein expressed in <i>E. coli</i> - lyophilised
Reconstitution	Reconstitute with 0.1 ml distilled water. For extended storage the addition of 0.1% BSA is recommended. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Buffer Solution	10mM Acetic Acid
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 0.1 ng/ug

Approx. Protein Concentrations	Total protein concentration 0.1 mg/ml after reconstitution
External Database Links	<p>UniProt: P13232 Related reagents</p> <p>Entrez Gene: 3574 IL7 Related reagents</p>
Product Information	<p>Recombinant human IL-7 is a 17.4 kDa protein comprised of 153 amino acids corresponding to the full-length mature sequence.</p> <p>IL-7 is a potent lymphoid cell growth factor acting on Pre-B, Pro-B and early T cells. In a dose dependent cell proliferation assay using 2E8 cells this product shows a specific activity of > 2 x 10⁶ units/mg.</p> <p>Recombinant human IL-7 is recognized by Rat anti Human IL-7 antibody clone BVD10-40F6 in ELISA.</p>
Protein Molecular Weight	17.5kDa (153 amino acids)
Purity	>98 % by SDS PAGE and HPLC
ELISA	Recombinant human IL-7 may be used as a standard for ELISA applications with either anti-Human IL-7 (purified) or anti-Human IL-7 (biotin conjugate) .
Western Blotting	Recombinant human IL-7 may be used as the positive control for western blotting applications with either anti-Human IL-7 (purified) or anti-Human IL-7 (biotin conjugate) .
Storage	<p>Prior to reconstitution store at -20°C. Following reconstitution store at -20°C.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PHP04610527
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376

Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

To

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

find a

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M362335:200501'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)