

Datasheet: PHP047A

BATCH NUMBER 155291

Description:	RECOMBINANT HUMAN INTERLEUKIN-10
Name:	IL-10
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			0.2 - 20ng/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate positive/negative controls.

Target Species	Human
Product Form	Purified recombinant IL-10 - lyophilized
Reconstitution	<p>Reconstitute with 0.1ml of 5mM sodium phosphate pH 7.2. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.</p> <p>After the initial reconstitution step, the product should be diluted to a working concentration using an aqueous buffer containing 0.1% BSA (or any albumin), or diluted in culture media containing serum such as FCS.</p>
Preparation	Recombinant protein expressed in <i>Escherichia coli</i>
Buffer Solution	Sodium Phosphate
Preservative Stabilisers	None present

Carrier Free	Yes
Endotoxin Level	< 0.1 ng/ug
Approx. Protein Concentrations	Total protein concentration 0.1mg/ml after reconstitution.
External Database Links	<p>UniProt: P22301 Related reagents</p> <p>Entrez Gene: 3586 IL10 Related reagents</p>
Product Information	<p>Recombinant Human interleukin-10 is a 161 amino acid recombinant IL-10 cytokine, with a range of reported properties including inhibition of macrophage mediated cytokine synthesis, suppression of the delayed type hypersensitivity response, and stimulation of the Th2 cells response.</p> <p>Recombinant human IL-10 has an ED₅₀ of <2.0ng/ml as determined by co-stimulation of MC/9 cells with human IL-4.</p>
Protein Molecular Weight	18.6 kD
Activity	5 x 10 ⁵ units/mg
Purity	>98% by SDS PAGE and HPLC analysis
Western Blotting	PHP047A may be used under reducing or non-reducing western blotting conditions with AHP795 .
Storage	<p>Prior to reconstitution store at -20°C. Following reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PHP047A 10527
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

RABBIT ANTI HUMAN INTERLEUKIN-10 (AHP795)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M414366:221206'

Printed on 25 Mar 2023

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)