

Datasheet: PHP048

BATCH NUMBER 169183

Description:	RECOMBINANT HUMAN INTERLEUKIN-11
Name:	IL-11
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA			▪	
Functional Assays	▪			0.1ng/ml - 10ng/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20ul distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Source	E.coli
Preservative Stabilisers	None present
Carrier Free	Yes
Approx. Protein Concentrations	0.5 mg/ml after reconstitution
External Database Links	UniProt:

[P20809](#) [Related reagents](#)

Entrez Gene:

[3589](#) IL11 [Related reagents](#)

Product Information **Human Interleukin 11** is a potent lymphoid cell growth factor which stimulates the growth and survival of certain B cells and T cells. This product has an ED50 of 0.5ng/ml as determined in a proliferation assay using murine T1165 hybridoma cells.

Protein Molecular Weight 19.5kDa (179 Amino acid sequence)

Purity >98% by SDS PAGE

Storage Prior to reconstitution store at -20°C. Following reconstitution store at -20°C. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10268 available at: <https://www.bio-rad-antibodies.com/SDS/PHP04810268>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com **Worldwide** Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com **Europe** Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M420081:230706'

Printed on 18 Nov 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)