

Datasheet: PHP147

BATCH NUMBER 160616

Description:	PURIFIED HUMAN SEX HORMONE BINDING GLOBULIN
Name:	SEX HORMONE BINDING GLOBULIN
Other names:	SHBG
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Functional Assays			▪	
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified human sex hormone binding globulin - liquid
Preparation	Purified human sex hormone binding globulin from human serum plasma.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Total protein concentration 1.0 mg/ml
External Database Links	UniProt: P04278 Related reagents

Entrez Gene:

[6462](#) SHBG [Related reagents](#)

Product Information **Purified Human sex hormone binding globulin**, is a ~90 kDa homodimeric glycoprotein composed of two 45 kDa subunits linked via disulphide bridges. It is an androgen transport protein which binds and regulates the free concentration of steroid hormones and is involved in steroid signalling at the cell membrane.

Protein Molecular Weight 45 kDa

Purity >90% by SDS PAGE

References

1. Rosner, W. *et al.* (1999) Sex hormone-binding globulin mediates steroid hormone signal transduction at the plasma membrane. [J Steroid Biochem Mol Biol. 69 \(1-6\): 481-5.](#)
2. Kahn, S.M. *et al.* (2002) Sex hormone-binding globulin is synthesized in target cells. [J Endocrinol. 175 \(1\): 113-20.](#)
3. Zheng, X. *et al.* (2015) Analysis of Hormone-Protein Binding in Solution by Ultrafast Affinity Extraction: Interactions of Testosterone with Human Serum Albumin and Sex Hormone Binding Globulin. [Anal Chem. 87 \(22\): 11187-94.](#)

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10209 available at:
10209: <https://www.bio-rad-antibodies.com/uploads/MSDS/10209.pdf>

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M369510:200529'

Printed on 08 Feb 2023
