

Datasheet: PHP149

BATCH NUMBER 159086

Description:	PURIFIED HUMAN ALBUMIN
Name:	SERUM ALBUMIN
Other names:	HSA
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	
Functional Assays			▪	
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species

Human

Product Form

Purified human serum albumin - lyophilized

Reconstitution

Reconstitute with phosphate buffer pH>7.0 containing 0.15M NaCl
 Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
 For long term storage the addition of 0.09% sodium azide is recommended.
 N.B. For functional studies do not add sodium azide

Preparation

Purified albumin from human urine

Preservative Stabilisers

None present

**External Database
Links**

UniProt:

[P02768](#) [Related reagents](#)

Entrez Gene:

[213](#) ALB [Related reagents](#)

Product Information **Purified human albumin** is a ~66 kDa protein that binds and transports steroids, fatty acids, metal ions and thyroid hormones. Albumin accounts for approximately 50% of blood serum protein.

Protein Molecular Weight 66 kDa

Purity >99% by SDS PAGE

References 1. Nicholson, J.P. *et al.* (2000) The role of albumin in critical illness. [Br J Anaesth. 85 \(4\): 599-610.](#)

Storage Prior to reconstitution store at +4°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10268 available at:
<https://www.bio-rad-antibodies.com/SDS/PHP14910268>

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M412118:221109'

Printed on 19 Jan 2024