

## Datasheet: PHP149

**BATCH NUMBER 171123**

<b>Description:</b>	PURIFIED HUMAN ALBUMIN
<b>Name:</b>	SERUM ALBUMIN
<b>Other names:</b>	HSA
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	1 mg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting			▪	
Functional Assays			▪	
Immunoassay	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Target Species

Human

#### Product Form

Purified human serum albumin - lyophilized

#### Reconstitution

Reconstitute with phosphate buffer pH>7.0 containing 0.15M NaCl

Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.

For long term storage the addition of 0.09% sodium azide is recommended.

N.B. For functional studies do not add sodium azide

#### Preparation

Purified albumin from human urine

#### Preservative Stabilisers

None present

**External Database  
Links**

**UniProt:**

[P02768](#)   [Related reagents](#)

**Entrez Gene:**

[213](#)   ALB   [Related reagents](#)

---

**Product Information**   **Purified human albumin** is a ~66 kDa protein that binds and transports steroids, fatty acids, metal ions and thyroid hormones. Albumin accounts for approximately 50% of blood serum protein.

---

**Protein Molecular Weight**   66 kDa

---

**Purity**   >99% by SDS PAGE

---

**References**   1. Nicholson, J.P. *et al.* (2000) The role of albumin in critical illness. [Br J Anaesth. 85 \(4\): 599-610.](#)

---

**Storage**   Prior to reconstitution store at +4°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee**   Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information**   Material Safety Datasheet documentation #10268 available at:  
<https://www.bio-rad-antibodies.com/SDS/PHP14910268>

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

---

**Regulatory**   For research purposes only

**North & South America**   Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M412118:221109'

**Printed on 25 Mar 2023**