

Datasheet: PHP161

| | |
|----------------------|-----------------------------------|
| Description: | RECOMBINANT HUMAN DEFENSIN BETA 2 |
| Name: | DEFENSIN BETA 2 |
| Other names: | BD-2 |
| Format: | Rec. Protein |
| Product Type: | Recombinant Protein |
| Quantity: | 20 µg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|-------------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | 0.2 - 0.4ng/well |
| Western Blotting | ▪ | | | 1.5 - 3.0ng/lane |
| Functional Assays | ▪ | | | 10 - 100ng/ml |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|---|
| Target Species | Human |
| Product Form | Purified recombinant protein - lyophilized |
| Reconstitution | Reconstitute with 0.2ml 10 mM acetic acid. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage the addition of 0.1% bovine serum albumin (BSA) is recommended. |
| Preparation | Purified recombinant BD-2, expressed in <i>E. coli</i> |
| Source | E.coli |
| Preservative Stabilisers | None present |
| Carrier Free | Yes |

| | |
|---------------------------------------|---|
| Endotoxin Level | < 0.1 ng/ug |
| Approx. Protein Concentrations | 0.1 mg/ml after reconstitution |
| External Database Links | <p>UniProt: O15263 Related reagents</p> <p>Entrez Gene: 1673 DEFB4A Related reagents</p> |
| Synonyms | DEFB102, DEFB2, DEFB4 |
| Product Information | <p>Recombinant human defensin beta-2 (BD-2), is a 4.3 kDa peptide comprised of 41 amino acid residues. BD-2, a member of the beta defensin family, is secreted at epithelial surfaces of the skin and respiratory tract and by some leukocytes. This defensin is induced by bacterial products and cytokines during inflammation and functions as part of the innate immune system, having a wide ranging antimicrobial activity. BD-2 also functions as a chemoattractant to immature dendritic cells and memory T cells and acts as a ligand for TLR4, upregulating co-stimulatory molecules and inducing dendritic cell maturation.</p> |
| Protein Molecular Weight | 4.3 kDa (41 amino acid residues) |
| Activity | Determined by its ability to chemoattract immature dendritic cells using a concentration of 10.0-100.0 ng/ml |
| Purity | >98% by SDS PAGE and HPLC |
| Amino Acid Sequence | GIGDPVTCLKSGAICHPVFCPRRYKQIGTCGLPGTKCKKP |
| ELISA | Recombinant human BD-2 may be used as a standard in ELISA applications with either a purified human BD-2 antibody (AHP849) or a biotinylated human BD-2 antibody (AHP489B). |
| Western Blotting | Recombinant human BD-2 may be used as the positive control in Western Blotting applications with either a purified human BD-2 antibody (AHP849) or a biotinylated human BD-2 antibody (AHP489B). |
| References | <ol style="list-style-type: none"> 1. Bals, R. <i>et al.</i> (1998) Human beta-Defensin 2 is a salt-sensitive peptide antibiotic expressed in human lung. J. Clin. Invest. 102: 874-880. 2. Yang, D. <i>et al.</i> (1999) Beta-defensins: linking innate and adaptive immunity through dendritic and T cell CCR6. Science 286: 525-528. |
| Storage | <p>Prior to reconstitution store at -20°C.</p> <p>After reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended.</p> |

This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

| | |
|------------------|---|
| Guarantee | Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date. |
|------------------|---|

| | |
|--------------------------------------|---|
| Health And Safety Information | Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PHP16110527 |
|--------------------------------------|---|

| | |
|-------------------|----------------------------|
| Regulatory | For research purposes only |
|-------------------|----------------------------|

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M433638:241118'

Printed on 18 Nov 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)