

Datasheet: PHP198

BATCH NUMBER 169994

Description:	RECOMBINANT HUMAN LEPTIN
Name:	LEPTIN
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 1.0 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage, the addition of 0.1% bovine serum albumin (BSA) is recommended.
Preparation	Purified recombinant leptin expressed in <i>E.coli</i> . Sterile filtered through a 0.2 micron filter
Source	E.coli
Buffer Solution	Lyophilized with 0.1% TFA
Preservative Stabilisers	None present

Carrier Free	Yes
Endotoxin Level	< 1.0 EU/ug
Approx. Protein Concentrations	1.0 mg/ml after reconstitution
External Database Links	<p>UniProt: P41159 Related reagents</p> <p>Entrez Gene: 3952 LEP Related reagents</p>
Synonyms	OB, OBS
Product Information	<p>Recombinant human leptin is a leptin preparation made in <i>E. coli</i>. Leptin is a major adipose-derived hormone and member of the leptin family, which plays a key role in appetite control and the regulation of energy expenditure.</p> <p>Secreted leptin associates with proteins within the bloodstream and is transported to the brain, where it inhibits or up-regulates the activity of neuropeptides such as neuropeptide Y (NPY) and alpha-melanocyte-stimulating hormone respectively. Leptin is a recognised marker of obesity, with a correlation between circulating levels of leptin and body fat mass, and defects in leptin production result in severe obesity. Regulators of leptin expression include insulin, cortisol, cAMP and thiazolidinediones.</p>
Protein Molecular Weight	16.0 kDa (147 amino acid residues)
Activity	PHP198 is biologically active in the <i>ob/ob</i> mouse obesity model. The <i>ob/ob</i> mice were treated via intraperitoneal injection once daily at a dose of 5 µg Leptin/gm of body weight for 7 days. Significant effects on body weight, food consumption, and plasma glucose levels were observed to saline-treated controls.
Purity	>98% by SDS PAGE and HPLC
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PHP198

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M433640:241118'

Printed on 29 Jan 2026
