

## Datasheet: PHP199

**BATCH NUMBER 156423**

<b>Description:</b>	RECOMBINANT HUMAN PAI-1
<b>Name:</b>	PAI-1
<b>Other names:</b>	SERPINE1
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	10 µg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 100 ul distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage, the addition of 0.1% BSA is recommended.
<b>Preparation</b>	Purified recombinant PAI-1 expressed in <i>E.coli</i> .
<b>Buffer Solution</b>	Sodium Acetate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Endotoxin Level</b>	< 1.0 EU/ug

<b>Approx. Protein Concentrations</b>	0.1 mg/ml after reconstitution
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P05121</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">5054</a>    SERPINE1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	PAI1, PLANH1
<b>Product Information</b>	PAI-1 (plasminogen activator inhibitor-1) is a serine protease inhibitor and member of the serpin family, which acts as the primary inhibitor of tissue plasminogen activator (tPA) and urokinase plasminogen activator (uPA), and hence fibrinolysis. Increased levels of PAI-1 are associated with diabetes, cardiovascular disease and an increased occurrence of thrombosis in obesity, and can be induced by obesity related factors including insulin, glucose and VLDL (very-low-density lipoprotein).
<b>Protein Molecular Weight</b>	42.7 kDa (379 amino acid residues).
<b>Activity</b>	Determined by its inhibitory effect against single chain tPA induced cleavage of a chromogenic substrate in Imidazole Buffer at 37°C. Half maximal inhibition against 1.0 µg/ml of single chain tPA was obtained at a concentration of 2.0 µg/ml.
<b>Purity</b>	>95% by SDS PAGE and HPLC analysis.
<b>Storage</b>	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Guarantee</b>	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10168 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PHP199">https://www.bio-rad-antibodies.com/SDS/PHP199</a> 10168
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M374911:201208'

