

## Datasheet: PHP201

**BATCH NUMBER 157263**

<b>Description:</b>	RECOMBINANT HUMAN VISFATIN
<b>Name:</b>	VISFATIN
<b>Other names:</b>	NAMPT
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	25 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 0.25 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. N.B. Due to solubility reasons, it is recommended that the solution be kept at low pH.
<b>Preparation</b>	Purified recombinant visfatin expressed in <i>E.coli</i>
<b>Buffer Solution</b>	10mM HCl
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Endotoxin Level</b>	< 1.0 EU/ug

<b>Approx. Protein Concentrations</b>	Total protein concentration 0.1 mg/ml after reconstitution
<b>External Database Links</b>	<b>UniProt:</b> <a href="#">P43490</a> <a href="#">Related reagents</a>  <b>Entrez Gene:</b> <a href="#">10135</a> NAMPT <a href="#">Related reagents</a>
<b>Synonyms</b>	PBEF, PBEF1
<b>Product Information</b>	<p>Visfatin is a putative cytokine originally identified as PBEF (Pre-B cell-enhancing factor), subsequently determined to be a nicotinamide phosphoribosyltransferase (NAMPT) enzyme and member of the NAPRTase family, involved in catalyzing the biosynthesis of NAD from nicotinamide.</p> <p>Visfatin is highly expressed in bone marrow, muscle and liver tissue and increases in expression in response to endotoxin, in amniotic epithelial cells following lipopolysaccharide (LPS) stimulation, and during the activation of T and B cells, dendritic cells and macrophages, suggesting a role in inflammatory processes.</p>
<b>Protein Molecular Weight</b>	52.5 kDa (465 amino acid residues)
<b>Activity</b>	The ED <sub>50</sub> was determined by the dose-dependant proliferation of the RPMI 8226 cells. The expected ED <sub>50</sub> for this effect is 15.0-20.0 ng/ml.
<b>Purity</b>	>98% by SDS PAGE and HPLC analysis.
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.</p> <p>After reconstitution store at -20°C.</p> <p>Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
<b>Health And Safety Information</b>	<p>Material Safety Datasheet documentation #10268 available at:</p> <p><a href="https://www.bio-rad-antibodies.com/SDS/PHP20110268">https://www.bio-rad-antibodies.com/SDS/PHP20110268</a></p>
<b>Regulatory</b>	For research purposes only

<b>North &amp; South America</b>	Tel: +1 800 265 7376	<b>Worldwide</b>	Tel: +44 (0)1865 852 700	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21
	Fax: +1 919 878 3751		Fax: +44 (0)1865 852 739		Fax: +49 (0) 89 8090 95 50
	Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>		Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>		Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M362371:200501'

**Printed on 19 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)