

Datasheet: PHP213

Description:	RECOMBINANT HUMAN BDNF
Name:	BDNF
Other names:	BRAIN-DERIVED NEUROTROPHIC FACTOR
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 100 ul distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified recombinant BDNF expressed in <i>E.coli</i>
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 1.0 EU/ug
Approx. Protein Concentrations	0.1 mg/ml after reconstitution

**External Database
Links**

UniProt:

[P23560](#) [Related reagents](#)

Entrez Gene:

[627](#) BDNF [Related reagents](#)

Product Information Brain-derived neurotrophic factor (BDNF) is a potent neurotrophin and member of the NGF-beta family, secreted from dense core vesicles, which promotes the growth, differentiation and survival of neurons within, or directly connected to, the central nervous system.

BDNF signals through binding to the high affinity receptor trkB, and also binds to the low affinity neurotrophin receptor LNGFR. BDNF is involved in the modulation of short-term synaptic function and in the formation of long-term memory, and reduced levels of BDNF are seen in patients with Alzheimer's and Huntington disease ([Girota et al. 2022](#)).

Protein Molecular Weight 27.0 kDa (119 amino acid homodimer)

Activity Determined by its ability to stimulate the proliferation of rat C6 cells. The expected ED₅₀ for this effect is 0.5-1.0 µg/ml.

Purity >98% by SDS PAGE and HPLC analysis

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.

Health And Safety Information Material Safety Datasheet documentation #10268 available at: 10268: <https://www.bio-rad-antibodies.com/uploads/MSDS/10268.pdf>

Regulatory For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M362622:200514'

Printed on 05 Apr 2022