

Datasheet: PHP236

BATCH NUMBER 162302

Description:	NATIVE HUMAN THYROGLOBULIN
Name:	THYROGLOBULIN
Other names:	Tg
Format:	Purified
Product Type:	Purified Protein
Quantity:	1 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human thyroid glands - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present.
Approx. Protein Concentrations	1.0 mg/ml
External Database Links	<p>UniProt: P01266 Related reagents</p> <p>Entrez Gene: 7038 TG Related reagents</p>

Product Information **Purified Human Thyroglobulin** is derived from human thyroid glands and is the precursor of thyroxine (T3) and triiodothyronine (T4). Thyroglobulin is a 2768 amino acid secreted glycoprotein with a 19 aa N-terminal signal peptide, made within and targeted to the thyroid.

Mutations of the TG gene may lead to the development of Thyroid dysmorphogenesis 3 ([TDH3](#)) a condition presenting with predominant goiters on the majority of patients and may lead to varying degrees of neurologic and intellectual defects. Variations in the TG gene may also result in increased susceptibility to the development of Autoimmune thyroid disease 3 ([AITD3](#)), an autoimmune disorder comprised of two related diseases of the thyroid: [Graves disease](#) and [Hashimoto thyroiditis](#) ([Caturegli et al. 2014](#)).

Protein Molecular Weight 660 kDa by chromatography and electrophoresis, under certain conditions the 330 kDa subunit may be observed.

Purity SDS PAGE: >98%, ultra-pure, IgG low

Storage Store at -20°C only.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10209 available at: <https://www.bio-rad-antibodies.com/SDS/PHP236>
10209

Donor material tested and found negative for HIV1 and 2 antibodies, HBsAg and HCV antibodies.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M395788:220519'

Printed on 16 May 2024