

Datasheet: PHP239

BATCH NUMBER 180403

Description:	NATIVE HUMAN LACTOFERRIN
Name:	LACTOFERRIN
Other names:	LACTOTRANSFERRIN
Format:	Purified
Product Type:	Purified Protein
Quantity:	5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using the appropriate negative/positive controls.

Target Species	Human
Product Form	Purified protein from human milk - lyophilised
Reconstitution	Reconstitute with 5 ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified by chromatography.
Preservative Stabilisers	None present

External Database Links

UniProt:

[P02788](https://www.uniprot.org/entry/P02788)

[Related reagents](#)

Entrez Gene:

[4057](https://www.ncbi.nlm.nih.gov/gene/4057)

LTF

[Related reagents](#)

Synonyms	LF
Product Information	<p>Human lactoferrin is a ~78 kDa globular protein found in body secretions such as milk and saliva and is a member of the transferrin family proteins. Transferrins control the level of free iron in the body, by binding and transporting iron. Lactoferrin does not appear to have a transporting function, instead it has antimicrobial activity, by keeping the levels of iron in body fluids low and thus preventing microbes from acquiring iron. Removing free iron from body fluids and inflamed areas also serve to protect against the damaging effects of free radicals.</p> <p>Lactoferrin is used as a therapeutic and prophylactic treatment of various illnesses, such as cancers and intestinal illnesses.</p>
Purity	~ 90% by Biuret method
References	<p>1. Wahl, A. <i>et al.</i> (2015) Breast Milk of HIV-Positive Mothers Has Potent and Species-Specific <i>In Vivo</i> HIV-Inhibitory Activity. J Virol. 89 (21): 10868-78.</p> <p>2. Yamashita, K. <i>et al.</i> (2017) Impact of a 7-Day Field Training on Oral Health Condition in Japan Ground Self-Defense Force Personnel. Mil Med. 182 (7): e1869-77.</p>
Storage	<p>Prior to reconstitution store at +4°C.</p> <p>After reconstitution store at +4°C.</p> <p>DO NOT FREEZE. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	12 months from date of despatch
Health And Safety Information	<p>Material Safety Datasheet documentation #10268 available at: https://www.bio-rad-antibodies.com/SDS/PHP23910268</p> <p>Donor material tested and found negative for HIV, HTLV, HCV and HBsAg.</p> <p>As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.</p>
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M371195:200602'

Printed on 18 Jan 2024