

Datasheet: PHP245

BATCH NUMBER 157776

Description:	RECOMBINANT HUMAN apo-SERUM AMYLOID A1
Name:	apo-SERUM AMYLOID A1
Other names:	apo-SAA1
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	50 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 0.5 ml distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. For extended storage, the addition of 0.1% BSA is recommended.
Preparation	Purified recombinant apo-SAA1 expressed in <i>E.coli</i>
Buffer Solution	5mM Tris, pH7.6
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 1.0 EU/ug

Approx. Protein Concentrations	0.1 mg/ml after reconstitution
External Database Links	<p>UniProt: P02735 Related reagents</p> <p>Entrez Gene: 6288 SAA1 Related reagents</p>
Product Information	Serum Amyloid A (apo-SAA) is an acute phase response protein and injury specific apolipoprotein, synthesized primarily in the liver, which is found in trace amounts in association with high density lipoproteins (HDL) in human serum, but increases dramatically during an acute response to inflammation, infection, trauma, and malignant neoplasia, to become the most abundant HDL apolipoprotein. The acute phase isoforms of apo-SAA are encoded by the genes SAA1 (Apo-SAA1) and SAA2 (Apo-SAA2), and both include allelic variants.
Protein Molecular Weight	11.7 kDa (105 amino acid residues)
Activity	Determined by its ability to chemoattract human monocytes using a range of 10.0-10.0 ng/ml.
Purity	>98% by SDS PAGE and HPLC analysis
Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PHP24510527
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M379728:210408'

Printed on 18 Jan 2024

