

## Datasheet: PHP279

**BATCH NUMBER 165069**

<b>Description:</b>	NATIVE HUMAN FACTOR H
<b>Name:</b>	FACTOR H
<b>Other names:</b>	COMPLEMENT FACTOR H
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Quantity:</b>	0.25 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Human
<b>Product Form</b>	Purified from human serum - liquid
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Current, batch-specific concentration 1.05 mg/ml
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P08603</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">3075</a>    CFH    <a href="#">Related reagents</a></p>

**Synonyms** HF, HF1, HF2

---

**Product Information** Human factor H is a single chain glycoprotein present in normal human serum at about 500 µg/ml that regulates the formation and function of complement C3 and C5 convertase enzymes. Regulatory activity is attributed to it's ability to recognize and bind C3b fragments.

---

**Protein Molecular Weight** 150 kDa

---

**Purity** >85%

---

**Storage** Store at -70°C.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information** Material Safety Datasheet documentation #10209 available at:  
<https://www.bio-rad-antibodies.com/SDS/PHP279>  
10209

The original starting material was tested and found to be negative for HBsAg, HCV and antibodies to HIV.

---

**Regulatory** For research purposes only

---

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M394734:220215'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)