

Datasheet: PHP283

BATCH NUMBER 165919

Description:	PURIFIED HUMAN C3
Name:	C3
Other names:	COMPLEMENT COMPONENT C3
Format:	Purified
Product Type:	Purified Protein
Quantity:	0.25 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified human C3 - liquid
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Current, batch-specific concentration 1.15 mg/ml
External Database Links	<p>UniProt:</p> <p>P01024 Related reagents</p> <p>Entrez Gene:</p> <p>718 C3 Related reagents</p>

Synonyms CPAMD1

Product Information C3 is a complement protein that is activated through cleavage by complement 3 (C3) convertases to C3a and C3b. These are further degraded into iC3b and C3dg/C3d. C3 is crucial in the induction of tolerance generated when an antigen is introduced into immunoprivileged sites and this is exploited by pathogens and cancer cells to evade the immune system by inhibiting complement activation.

Purity >85% by SDS PAGE

Storage Store at -70°C.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee Guaranteed until date of expiry. Please see product label.

Health And Safety Information Material Safety Datasheet documentation #10209 available at:
<https://www.bio-rad-antibodies.com/SDS/PHP283>
10209

Donor material tested and found negative for HIV1, HIV2, HBsAg and HCV.

As no test can completely guarantee this material to be free of pathogens it should be handled as potentially infectious.

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M395837:220525'

Printed on 12 Feb 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)