

Datasheet: PHP293

BATCH NUMBER 153177

Description:	RECOMBINANT HUMAN VEGF-A
Name:	VEGF-A
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	0.2 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Human
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Centrifuge vial prior to reconstitution. Reconstitute to 0.5 mg/ml by adding 0.4 ml ddH ₂ O. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Mix gently after reconstitution. Do not vortex.
Preparation	Recombinant protein expressed in <i>E. coli</i> and purified by affinity chromatography
Buffer Solution	20 mM Phosphate Buffer, 0.1 M Sodium Chloride
Preservative Stabilisers	1% Trehalose
Approx. Protein Concentrations	0.5 mg/ml after reconstitution
External Database Links	UniProt:

Product Information	Recombinant Human VEGF-A <p>Vascular endothelial growth factor alpha forms a homodimer and is a potent stimulator of angiogenesis of both normal and cancerous cells. It acts as a regulator of vasculogenesis, angiogenesis and endothelial cell growth. VEGF-A is secreted by many different cell types such as endothelial cells (Nissen <i>et al.</i> 1998), smooth muscle cells (Brogi <i>et al.</i> 1994), neutrophils (Gaudry <i>et al.</i> 1997), platelets (Banks <i>et al.</i> 1998), macrophages and 60% of all tumours (Berse <i>et al.</i> 1992).</p> <p>The sequence of this recombinant protein product is amino acid (aa) 27 - aa 141 + CDKPRR but deficient from aa 142-226, and it is closest to isoform 121. Isoform 121 and 111 lack exons 6 and 7, and for this reason neither of these isoforms bind to the extracellular matrix (Krillicke <i>et al.</i> 2017). This product has been demonstrated for use, in a pharmacokinetic (PK) antigen capture ELISA with the monoclonal antibody drugs ranibizumab and bevacizumab using anti-ranibizumab antibodies for detection.</p>
Protein Molecular Weight	Predicted 14 kDa. The apparent molecular mass of the human VEGF-A monomer is approximately 17 kDa determined by SDS-PAGE under reducing conditions, and 17 kDa (monomer) and 30 kDa (dimer) under non-reducing conditions.
Purity	≥90% determined by SDS-PAGE under reducing conditions and visualized by coomassie blue staining
Amino Acid Sequence	APMAEGGGQNHHEVVKFMDVYQRSYCHPIETLVDIFQEYYPDEIEYIFKPSCVPLMRCGGCCNDEGLE CVPTEESNITMQIMRIKPHQGQHIGEMSFLQHNKCECRPKKDRARQEKC DKPRR
Storage	Prior to reconstitution store at -70°C. Following reconstitution store at -20°C. This product should be stored undiluted. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein.
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #20395 available at: https://www.bio-rad-antibodies.com/SDS/PHP29320395
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[HUMAN ANTI RANIBIZUMAB \(DRUG/TARGET COMPLEX\) \(HCA304\)](#)

[HUMAN ANTI RANIBIZUMAB \(HCA307\)](#)

North & South Tel: +1 800 265 7376

Worldwide Tel: +44 (0)1865 852 700

Europe Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M350048:190307'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)