

Datasheet: PIP005

BATCH NUMBER 163150

Description:	NATIVE CYTOMEGALOVIRUS
Name:	CYTOMEGALOVIRUS
Other names:	CMV
Format:	Inactivated Pathogen
Product Type:	Antigen
Quantity:	1 ml

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Viral
Product Form	Inactivated Cytomegalovirus - liquid
Preparation	Cytomegalovirus strain AD-169, cultured in MRC-5, a diploid cell line of human lung origin. MRC-5 monolayers are infected in a manner which will ensure that the full host of viral proteins are present upon harvest. Optimally infected tissue is collected in glycine buffered saline. The cells are disrupted to release cellular CMV antigen and cellular debris is removed by centrifugation. The antigen preparation is inactivated using gamma radiation, which primarily damages viral genetic material.
Buffer Solution	Glycine buffered saline
Preservative Stabilisers	None present
Approx. Protein Concentrations	Current, batch-specific concentration 3.5 mg/ml

Product Information Cytomegalovirus is a member of the *Herpesviridae* family. It may be transmitted from person to person through urine, saliva, other bodily fluids, breast milk and transplanted organs. It is a common infection which mostly passes unnoticed. The infection may be dangerous to people who are immunocompromised such as those who are HIV positive or organ transplant recipients. The infection may also be dangerous to pre- and post-natal infants. Antiviral drugs are useful in the treatment of CMV infection.

Activity Antigenic activity is 129% of internal reference standard.

Instructions For Use PIP005 should be sonicated immediately before use to ensure the preparation is uniform. The product may be used in a variety of immunoassay formats or may be further purified to meet the requirements of a particular assay format.

Storage Store at -70°C.
Storage in frost-free freezers is not recommended.
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10220 available at: <https://www.bio-rad-antibodies.com/SDS/PIP005>
10220
This product has been rendered inactive by standard procedures. However this material should still be handled as infectious and should be disposed of appropriately.

Regulatory For research purposes only

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M390480:210913'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)