

## Datasheet: PIP005

**BATCH NUMBER 163452**

<b>Description:</b>	NATIVE CYTOMEGALOVIRUS
<b>Name:</b>	CYTOMEGALOVIRUS
<b>Other names:</b>	CMV
<b>Format:</b>	Inactivated Pathogen
<b>Product Type:</b>	Antigen
<b>Quantity:</b>	1 ml

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	■			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Target Species

Viral

#### Product Form

Inactivated Cytomegalovirus - liquid

#### Preparation

Cytomegalovirus strain AD-169, cultured in MRC-5, a diploid cell line of human lung origin. MRC-5 monolayers are infected in a manner which will ensure that the full host of viral proteins are present upon harvest. Optimally infected tissue is collected in glycine buffered saline. The cells are disrupted to release cellular CMV antigen and cellular debris is removed by centrifugation. The antigen preparation is inactivated using gamma radiation, which primarily damages viral genetic material.

#### Buffer Solution

Glycine buffered saline

#### Preservative Stabilisers

None present

#### Approx. Protein Concentrations

Current, batch-specific concentration 3.3 mg/ml

**Product Information** Cytomegalovirus is a member of the *Herpesviridae* family. It may be transmitted from person to person through urine, saliva, other bodily fluids, breast milk and transplanted organs. It is a common infection which mostly passes unnoticed. The infection may be dangerous to people who are immunocompromised such as those who are HIV positive or organ transplant recipients. The infection may also be dangerous to pre- and post-natal infants. Antiviral drugs are useful in the treatment of CMV infection.

---

**Activity** Antigenic activity is 112% of internal reference standard.

---

**Instructions For Use** PIP005 should be sonicated immediately before use to ensure the preparation is uniform. The product may be used in a variety of immunoassay formats or may be further purified to meet the requirements of a particular assay format.

---

**Storage** Store at -70°C.  
Storage in frost-free freezers is not recommended.  
This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

---

**Guarantee** Guaranteed until date of expiry. Please see product label.

---

**Health And Safety Information** Material Safety Datasheet documentation #10220 available at: <https://www.bio-rad-antibodies.com/SDS/PIP005>  
10220  
This product has been rendered inactive by standard procedures. However this material should still be handled as infectious and should be disposed of appropriately.

---

**Regulatory** For research purposes only

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M407144:221004'

**Printed on 18 Jan 2024**

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)