

## Datasheet: PIP021

**BATCH NUMBER 160441**

<b>Description:</b>	NATIVE INFLUENZA A H1N1 ANTIGEN
<b>Name:</b>	INFLUENZA A H1N1 ANTIGEN
<b>Format:</b>	Inactivated Pathogen
<b>Product Type:</b>	Antigen
<b>Quantity:</b>	1 ml

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Haemagglutination	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Viral
<b>Product Form</b>	Inactivated Influenza A H1N1 antigen - liquid
<b>Preparation</b>	Influenza strain A/New Caledonia/20/99 (H1N1) cultured in embryonated chicken eggs. Eggs are infected by injection into the allantoic chamber, incubated for 72 hours, and refrigerated overnight prior to harvest. Allantoic fluids are harvested from the eggs and pooled. The antigen preparation is inactivated by gamma radiation.
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Current, batch-specific concentration 2.1 mg/ml
<b>Product Information</b>	<b>Native influenza H1N1 antigen preparation</b> is presented as semi-purified allantoic fluid from chicken eggs infected with New Caladonia 20/99 (H1N1) influenza A virus. Influenza

is an acute infection of the respiratory tract caused by the Influenza virus, of which there are several different types. Influenza A causes flu in birds and some mammals and is the most virulent of the influenza types, resulting in the most severe disease. Influenza infection in humans causes chills, fever, pharyngitis, muscle pain, headaches, coughing and weakness. The virus is transmitted through the air by coughs or sneezes or through saliva, nasal secretions, faeces and blood. Infection may occur through direct contact with infected bodily fluid or contaminated surfaces. Influenza viruses mutate rapidly, making vaccines less effective from one season to the next.

This particular strain of influenza is a reassortant virus that includes genetic components of the H1N1 influenza pandemic swine flu in a culture-adapted virus that has been selected for its qualities as a vaccine-producing strain. The live reassortant virus has passed safety testing in ferrets, in which it is attenuated for pathogenicity relative to wild-type A/California/4/2009 virus. Because A/California/4/2009 is identified as a variant strain of a pandemic virus, live derivatives must be handled under BSL2<sup>+</sup> facilities and procedures.

Native influenza H1N1 antigen preparation (PIP021) contains a high concentration of viral antigens as well as some egg proteins from allantoic fluid.

Hemagglutination endpoint assay; 30720 HA Units/ml

---

<b>Activity</b>	40960 HA units/ml as determined by hemagglutination endpoint assay.
-----------------	---

---

<b>Instructions For Use</b>	A precipitate may form upon thawing, due to high product concentration. Dilute prior to clarification or other manipulation.
-----------------------------	--

---

<b>References</b>	1. Sakamoto, H. <i>et al.</i> (2015) A Novel Optical Biosensing System Using Mach-Zehnder-Type Optical Waveguide for Influenza Virus Detection. <a href="#">Appl Biochem Biotechnol. Oct 24. [Epub ahead of print]</a>
-------------------	--

---

<b>Storage</b>	Store at -70°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.
----------------	---

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
------------------	--

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10518 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PIP021">https://www.bio-rad-antibodies.com/SDS/PIP021</a> 10518  This product has been rendered inactive by gamma irradiation. However this material should still be handled as infectious and should be disposed of appropriately.
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**North & South** Tel: +1 800 265 7376

**Worldwide**

Tel: +44 (0)1865 852 700

**Europe**

Tel: +49 (0) 89 8090 95 21

To

**America** Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

find a

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M393985:220125'

**Printed on 13 Aug 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)