

## Datasheet: PIP050A

**BATCH NUMBER 149345**

<b>Description:</b>	RECOMBINANT DENGUE VIRUS TYPE 4 NS1 ANTIGEN
<b>Name:</b>	DENGUE VIRUS TYPE 4 NS1 ANTIGEN
<b>Other names:</b>	DENV4 NS1 ANTIGEN
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	100 µg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Viral
<b>Product Form</b>	Purified recombinant protein with a C-terminal His-tag - liquid
<b>Preparation</b>	Recombinant dengue virus serotype 4 NS1 protein, sequence strain Dominica/814669 /1981, expressed in 293 human cells
<b>Buffer Solution</b>	Tris buffered saline
<b>Preservative Stabilisers</b>	None present
<b>Approx. Protein Concentrations</b>	Total protein concentration 0.726 mg/ml

#### Specificity

**Recombinant dengue virus type 4 NS1 antigen** is a purified preparation of the dengue virus serotype 4 (DENV4), non-structural protein 1 (NS1). DENV4 is one of four (DENV1-4) antigenically distinct, closely related viral serotypes belonging to the *Flaviviridae* family, genus flavivirus. All four serotypes are known to cause viral infection in

humans.

Dengue virus is common mosquito-borne infection and a major cause of morbidity in tropical and subtropical regions. There is currently no vaccine to prevent, or effective anti-viral drugs to treat, dengue virus infection. In many cases infection is asymptomatic and the majority of individuals who get ill only suffer the mild, non-specific febrile symptoms characteristic of dengue fever (DF). Only a minority of infections result in severe disease, manifesting as dengue hemorrhagic fever (DHF) or dengue shock syndrome (DSS). Dengue virus infection gives lifelong immunity to the serotype in question but subsequent infection with another serotype may increase the likelihood of severe disease.

The NS1 glycoprotein is essential for viral replication and viability, and since this protein is secreted into the bloodstream, tests have been developed to diagnose DENV infections using NS1, including antigen-capture ELISA, lateral flow antigen detection, and the measurement of NS1-specific IgM and IgG responses ([Guzman, M.G. et al. 2010](#)).

Recombinant DENV4 NS1 antigen is presented in its native folded state complete with post-translational modifications, delivering optimal antigenicity and making it suitable for use in vaccine research and serology-based assays.

---

<b>Purity</b>	>90% by SDS PAGE analysis
---------------	---------------------------

---

<b>References</b>	1. Rose, N. <i>et al.</i> (2016) Generation of monoclonal antibodies against native viral proteins using antigen-expressing mammalian cells for mouse immunization. <a href="#">BMC Biotechnol. 16 (1): 83.</a>
-------------------	---

---

<b>Further Reading</b>	1. Guzman, M.G. <i>et al.</i> (2010) Dengue: a continuing global threat. <a href="#">Nat Rev Microbiol. 8 (12 Suppl): S7-16.</a>
------------------------	--

---

<b>Storage</b>	Store at -70°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Acknowledgements</b>	His-tag is a registered trademark of EMD Biosciences.
-------------------------	---

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10560 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PIP050A">https://www.bio-rad-antibodies.com/SDS/PIP050A</a> 10560
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

**North & South America** Tel: +1 800 265 7376  
Fax: +1 919 878 3751

**Worldwide** Tel: +44 (0)1865 852 700  
Fax: +44 (0)1865 852 739

**Europe** Tel: +49 (0) 89 8090 95 21  
Fax: +49 (0) 89 8090 95 50

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M378973:210319'

**Printed on 05 Mar 2024**

---