

Datasheet: PMP24

BATCH NUMBER 165716

Description:	RECOMBINANT MOUSE TNF ALPHA
Name:	TNF ALPHA
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	20 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			0.05ng/ml - 20ng/ml

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified protein produced in <i>E. coli</i> - lyophilized
Reconstitution	Reconstitute with 0.2 ml distilled water. For extended storage the addition of 0.1% BSA is recommended. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Source	E.coli
Buffer Solution	10mM Sodium Phosphate
Preservative Stabilisers	None present
Carrier Free	Yes

Endotoxin Level	< 0.1 ng/ug
Approx. Protein Concentrations	0.1 mg/ml after reconstitution
External Database Links	<p>UniProt: P06804 Related reagents</p> <p>Entrez Gene: 21926 Tnf Related reagents</p>
Synonyms	Tnfa, Tnfsf2
Product Information	Tumor necrosis factor alpha (TNFα) is a potent lymphoid factor which exerts cytotoxic effects on a wide range of tumor cells and certain other target cells.
Protein Molecular Weight	17.3 kDa (156 amino acid residues)
Activity	1 x 10 ⁷ units/mg
Purity	>98% by SDS PAGE and HPLC analysis
ELISA	Recombinant mouse TNFα may be used as the standard for ELISA applications with anti-Mouse TNFα (purified) (AAM19GA).
Western Blotting	Recombinant mouse TNFα may be used as the positive control for western blotting applications with anti-Mouse TNFα (AAM19GA) .
References	<ol style="list-style-type: none"> 1. Dunzendorfer, S. <i>et al.</i> (2004) Toll-like receptor 4 functions intracellularly in human coronary artery endothelial cells: roles of LBP and sCD14 in mediating LPS responses. FASEB J. 2004 Jul;18(10):1117-9. 2. Rosowski, E.E. <i>et al.</i> (2011) Strain-specific activation of the NF-κB pathway by GRA15, a novel <i>Toxoplasma gondii</i> dense granule protein. J Exp Med. 208: 195-212.
Storage	<p>Store at -20°C only.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	<p>Material Safety Datasheet documentation #10527 available at: https://www.bio-rad-antibodies.com/SDS/PMP24 10527</p>

North & South Tel: +1 800 265 7376

Worldwide

Tel: +44 (0)1865 852 700

Europe

Tel: +49 (0) 89 8090 95 21

America Fax: +1 919 878 3751

Fax: +44 (0)1865 852 739

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_us@bio-rad.com

Email: antibody_sales_uk@bio-rad.com

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M415788:230113'

Printed on 18 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)