

Datasheet: PMP35

BATCH NUMBER 160458

Description:	RECOMBINANT MOUSE MCP-1
Name:	MCP-1
Other names:	CCL2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			1.0 - 20ng/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the protein for use in their own system using appropriate positive/negative controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 0.1ml distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein prepared from <i>E. coli</i> .
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 0.1 ng/ug

Approx. Protein Concentrations	0.1mg/ml after reconstitution
External Database Links	<p>UniProt: P10148 Related reagents</p> <p>Entrez Gene: 20296 Ccl2 Related reagents</p>
Synonyms	Je, Mcp1, Scya2
Product Information	Recombinant mouse MCP-1 is 13.8 kDa peptide corresponding to full-length mature mouse monocyte chemotactic protein 1 . MCP-1, otherwise known as CCL2 or JE, is a secreted chemokine involved in the recruitment of monocytes to sites of injury or inflammation (Lu, B. et al., 1998).
Protein Molecular Weight	13.8kD (125 amino acid sequence)
Purity	>98% by SDS PAGE and HPLC analysis
ELISA	This product may be used as a standard for ELISA applications with our purified mouse MCP-1 antibody (AAM43).
References	1. Ren, Z. <i>et al.</i> (2019) Minocycline attenuates ethanol-induced cell death and microglial activation in the developing spinal cord. Alcohol. 79: 25-35.
Storage	<p>Prior to reconstitution store at -20°C. Following reconstitution store at -20°C.</p> <p>This product should be stored undiluted</p> <p>Storage in frost free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10268 available at: https://www.bio-rad-antibodies.com/SDS/PMP35 10268
Regulatory	For research purposes only

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M391194:211008'

Printed on 18 Nov 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)