

## Datasheet: PMP52

**BATCH NUMBER 147622**

<b>Description:</b>	MOUSE IgM KAPPA
<b>Name:</b>	IgM KAPPA
<b>Format:</b>	Purified
<b>Product Type:</b>	Purified Protein
<b>Clone:</b>	MADNP-5
<b>Isotype:</b>	IgM
<b>Quantity:</b>	1 mg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			1.25ug/ml

Where this protein has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the protein for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified monoclonal IgM - liquid
<b>Preparation</b>	Purified IgM prepared by affinity chromatography from tissue culture supernatant
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	0.1% Sodium Azide
<b>Approx. Protein Concentrations</b>	IgM concentration 1.0 mg/ml

### External Database Links

#### UniProt:

[P01837](#)     [Related reagents](#)  
[P01872](#)     [Related reagents](#)

**Entrez Gene:**

[16071](#) Igk-C [Related reagents](#)

[16019](#) Ighm [Related reagents](#)

[16019](#) Ighm [Related reagents](#)

---

**RRID** AB\_324175

---

**Purity** Single band by agarose electrophoresis and immunoelectrophoresis.

---

**References**

1. Hayer, S. *et al.* (2008) B-cell infiltrates induce endosteal bone formation in inflammatory arthritis. [J Bone Miner Res. 23: 1650-60.](#)
2. Rohrer, J.W. *et al.* (2006) Identification of oncofetal antigen/immature laminin receptor protein epitopes that activate BALB/c mouse OFA/iLRP-specific effector and regulatory T cell clones. [J Immunol. 176: 2844-56.](#)

---

**Storage** Store at +4°C or at -20°C if preferred.

This product should be stored undiluted.

Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.

---

**Guarantee** 12 months from date of despatch

---

**Health And Safety Information** Material Safety Datasheet documentation #10335 available at: <https://www.bio-rad-antibodies.com/SDS/PMP52>  
10335

---

**Regulatory** For research purposes only

---

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M369573:200529'

Printed on 12 Aug 2024