

Datasheet: PMP80

BATCH NUMBER 169866

Description:	RECOMBINANT MOUSE M-CSF
Name:	M-CSF
Other names:	MACROPHAGE COLONY-STIMULATING FACTOR
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	10 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 100 ul distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified recombinant M-CSF expressed in <i>E.coli</i>
Source	E.coli
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 1.0 EU/ug

Approx. Protein Concentrations	0.1 mg/ml after reconstitution
External Database Links	<p>UniProt: P07141 Related reagents</p> <p>Entrez Gene: 12977 Csf1 Related reagents</p>
Synonyms	Csfm
Product Information	<p>M-CSF (macrophage colony-stimulating factor) is a potent lineage restricted hematopoietic growth factor and member of the CSF family of hormone-like glycoproteins, which plays a key role in the regulation of cellular differentiation, proliferation, and the survival of blood monocytes, tissue macrophages and their progenitor cells.</p> <p>Along with sRANKL/CD254/OPGL, M-CSF is essential for osteoclastogenesis, acting through the tyrosine kinase receptor c-fms (colony stimulating factor 1 receptor). A genetic defect in M-CSF is responsible for osteopetrosis in mice (op/op), which are severely deficient in mature macrophages and have a restricted capacity for bone remodeling.</p>
Protein Molecular Weight	36.4 kDa (156 amino acid homodimer)
Activity	The ED ₅₀ determined by the dose-dependent stimulation of the proliferation of murine M-NFS-60 cells is < 1.0 ng/ml, corresponding to a specific activity of > 1 x 10 ⁶ units/mg.
Purity	≥ 98% by SDS PAGE and HPLC analysis
Storage	<p>Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
Health And Safety Information	Material Safety Datasheet documentation #10268 available at: https://www.bio-rad-antibodies.com/SDS/PMP80
Regulatory	For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M433646:241118'

