

## Datasheet: PMP80

**BATCH NUMBER 171138**

<b>Description:</b>	RECOMBINANT MOUSE M-CSF
<b>Name:</b>	M-CSF
<b>Other names:</b>	MACROPHAGE COLONY-STIMULATING FACTOR
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	10 µg

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Mouse
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 100 ul distilled water Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Purified recombinant M-CSF expressed in <i>E.coli</i>
<b>Source</b>	E.coli
<b>Preservative Stabilisers</b>	None present
<b>Carrier Free</b>	Yes
<b>Endotoxin Level</b>	< 1.0 EU/ug

<b>Approx. Protein Concentrations</b>	0.1 mg/ml after reconstitution
<b>External Database Links</b>	<p><b>UniProt:</b>  <a href="#">P07141</a>    <a href="#">Related reagents</a></p> <p><b>Entrez Gene:</b>  <a href="#">12977</a> Csf1    <a href="#">Related reagents</a></p>
<b>Synonyms</b>	Csfm
<b>Product Information</b>	<p>M-CSF (macrophage colony-stimulating factor) is a potent lineage restricted hematopoietic growth factor and member of the CSF family of hormone-like glycoproteins, which plays a key role in the regulation of cellular differentiation, proliferation, and the survival of blood monocytes, tissue macrophages and their progenitor cells.</p> <p>Along with sRANKL/CD254/OPGL, M-CSF is essential for osteoclastogenesis, acting through the tyrosine kinase receptor c-fms (colony stimulating factor 1 receptor). A genetic defect in M-CSF is responsible for osteopetrosis in mice (op/op), which are severely deficient in mature macrophages and have a restricted capacity for bone remodeling.</p>
<b>Protein Molecular Weight</b>	36.4 kDa (156 amino acid homodimer)
<b>Activity</b>	The ED <sub>50</sub> determined by the dose-dependent stimulation of the proliferation of murine M-NFS-60 cells is < 1.0 ng/ml, corresponding to a specific activity of > 1 x 10 <sup>6</sup> units/mg.
<b>Purity</b>	≥ 98% by SDS PAGE and HPLC analysis
<b>Storage</b>	<p>Prior to reconstitution store at -20°C.  After reconstitution store at -20°C.  Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
<b>Guarantee</b>	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10268 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PMP80">https://www.bio-rad-antibodies.com/SDS/PMP80</a>
<b>Regulatory</b>	For research purposes only

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M433646:241118'

