

Datasheet: PMP81

BATCH NUMBER 166834

Description:	RECOMBINANT MOUSE INTERLEUKIN-4
Name:	IL-4
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	20 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			0.2 - 0.4ng/well
Western Blotting	▪			1.5 - 3.0ng/lane
Functional Assays	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Mouse
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 200 ul distilled water. Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Purified recombinant IL-4 expressed in <i>E.coli</i>
Source	E.coli
Preservative Stabilisers	None present
Carrier Free	Yes
Endotoxin Level	< 1.0 EU/ug

Approx. Protein Concentrations	0.1 mg/ml after reconstitution
External Database Links	<p>UniProt: P07750 Related reagents</p> <p>Entrez Gene: 16189 IL4 Related reagents</p>
Synonyms	IL-4
Product Information	Recombinant Murine interleukin-4 is an <i>E.coli</i> derived preparation of Interleukin-4 (IL-4), also known as B-cell stimulatory factor 1. Interleukin-4 is a pleiotropic cytokine which interacts with multiple cell lineages, including T cells, fibroblasts and haematopoietic cells, and is one of the most important factors involved in B cell IgG1 and IgE class switching. IL-4 acts as a regulator of low affinity Fc receptor for IgE (CD23) expression on both lymphocytes and monocytes.
Protein Molecular Weight	13.5 kDa (120 amino acid residues)
Activity	The ED ₅₀ was determined by the dose-dependent proliferation of murine HT-2 cells is < 2.0 ng/ml, corresponding to a specific activity of > 5 x 10 ⁵ units/mg.
Purity	>98% by SDS PAGE and HPLC
ELISA	PMP81 may be used as a standard in an ELISA application with AAM36 .
Western Blotting	PMP81 may be used as a positive control for western blotting applications with AAM36 .
References	<ol style="list-style-type: none"> 1. Thomas, M.J. <i>et al.</i> (2002) CD8 T cells inhibit IgE via dendritic cell IL-12 induction that promotes Th1 T cell counter-regulation. J Immunol. 168 (1): 216-23. 2. Bénateau, M. <i>et al.</i> (2012) Combination of glycolysis inhibition with chemotherapy results in an antitumor immune response. Proc Natl Acad Sci U S A. 109: 20071-6. 3. Lam, R.S. <i>et al.</i> (2016) Unprimed, M1 and M2 Macrophages Differentially Interact with <i>Porphyromonas gingivalis</i>. PLoS One. 11 (7): e0158629.
Storage	<p>Prior to reconstitution store at -20°C. Following reconstitution store at -20°C.</p> <p>This product should be stored undiluted.</p> <p>Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.</p>
Guarantee	Guaranteed for 3 months from the date of reconstitution or until the date of expiry, whichever comes first. Please see label for expiry date.

**Health And Safety
Information**

Material Safety Datasheet documentation #10268 available at:
<https://www.bio-rad-antibodies.com/SDS/PMP81>
10268

Regulatory

For research purposes only

**North & South
America** Tel: +1 800 265 7376
Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M420161:230706'

Printed on 19 Jan 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)