

Datasheet: PPP022

**BATCH NUMBER 168289**

<b>Description:</b>	RECOMBINANT PIG INTERFERON GAMMA
<b>Name:</b>	IFN GAMMA
<b>Other names:</b>	IFNG
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Source</b>	<i>Pichia pastoris</i>
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose

**External Database  
Links**

**UniProt:**

[P17803](#)    [Related reagents](#)

**Entrez Gene:**

[396991](#) IFNG    [Related reagents](#)

---

**Product Information**    **Recombinant Pig interferon gamma**, (IFN $\gamma$ ), known previously as macrophage activating factor, belongs to the type II interferon family and is produced by T lymphocytes and NK cells.

Native pig interferon gamma is a 166 amino acid cytokine with potent activating, anti-viral and anti-proliferative properties.

IFN $\gamma$  is predominantly secreted by activated T lymphocytes in response to specific mitogens as a result of infection and has numerous immunological regulatory functions. Properties include macrophage and NK cell activation, upregulation of MHC class I and II molecules on various cells and induction of cytotoxicity.

---

**Protein Molecular Weight**    16.8 kDa

---

**Purity**    >95% by SDS PAGE analysis

---

**ELISA**    This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2390](#) as a coating antibody and [AHP2390B](#) as a detection antibody.

---

**References**    1. Dawson, H.D. *et al.* (2020) Porcine cytokines, chemokines and growth factors: 2019 update. [Res Vet Sci. 131: 266-300.](#)

---

**Storage**    Prior to reconstitution store at -20°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

---

**Guarantee**    12 months from date of despatch

---

**Health And Safety Information**    Material Safety Datasheet documentation #10591 available at:  
<https://www.bio-rad-antibodies.com/SDS/PPP02210591>

---

**Regulatory**    For research purposes only

---

**North & South America**    Tel: +1 800 265 7376

Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M415804:230113'

