

Datasheet: PPP025

BATCH NUMBER 156776

| | |
|----------------------|--|
| Description: | RECOMBINANT PIG INTERLEUKIN-10 |
| Name: | IL-10 |
| Other names: | CYTOKINE STIMULATING INHIBITORY FACTOR |
| Format: | Rec. Protein |
| Product Type: | Recombinant Protein |
| Quantity: | 5 µg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Western Blotting | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Pig |
| Product Form | Purified recombinant protein - lyophilized |
| Reconstitution | Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. |
| Preparation | Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 10% Trehalose |
| External Database Links | UniProt: |

Entrez Gene:

[397106](#) IL10 [Related reagents](#)

| | |
|--------------------------------------|---|
| Product Information | Recombinant Pig interleukin-10 (IL-10) , also known as cytokine synthesis inhibitory factor (CSIF), is a secreted anti inflammatory cytokine produced primarily by monocytes. This product is a purified recombinant cytokine corresponding to aa 18-175 of porcine IL-10. IL-10 inhibits the synthesis of pro-inflammatory cytokines such as IFN γ , IL-2, IL-4, GM-CSF and TNF α . |
| Protein Molecular Weight | 18.1 kDa |
| Purity | >95% by SDS PAGE analysis |
| ELISA | This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2393 as a coating antibody and AHP2393B as a detection antibody. |
| Storage | Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein. |
| Guarantee | 12 months from date of despatch |
| Health And Safety Information | Material Safety Datasheet documentation #10591 available at: https://www.bio-rad-antibodies.com/SDS/PPP025 10591 |
| Regulatory | For research purposes only |

North & South America Tel: +1 800 265 7376
Fax: +1 919 878 3751
Email: antibody_sales_us@bio-rad.com

Worldwide Tel: +44 (0)1865 852 700
Fax: +44 (0)1865 852 739
Email: antibody_sales_uk@bio-rad.com

Europe Tel: +49 (0) 89 8090 95 21
Fax: +49 (0) 89 8090 95 50
Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M334929:181203'

Printed on 22 Feb 2024