

Datasheet: PPP025

Description:	RECOMBINANT PIG INTERLEUKIN-10
Name:	IL-10
Other names:	CYTOKINE STIMULATING INHIBITORY FACTOR
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Source	<i>Pichia pastoris</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

**External Database
Links**

UniProt:

[Q29055](#) [Related reagents](#)

Entrez Gene:

[397106](#) IL10 [Related reagents](#)

Product Information **Recombinant Pig interleukin-10 (IL-10)**, also known as cytokine synthesis inhibitory factor (CSIF), is a secreted anti inflammatory cytokine produced primarily by monocytes. This product is a purified recombinant cytokine corresponding to aa 18-175 of porcine IL-10. IL-10 inhibits the synthesis of pro-inflammatory cytokines such as IFN γ , IL-2, IL-4, GM-CSF and TNF α .

Protein Molecular Weight 18.1 kDa

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2393](#) as a coating antibody and [AHP2393B](#) as a detection antibody.

References 1. Dawson, H.D. *et al.* (2020) Porcine cytokines, chemokines and growth factors: 2019 update. [Res Vet Sci. 131: 266-300.](#)

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at: <https://www.bio-rad-antibodies.com/SDS/PPP025>
10591

Regulatory For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M415807:230113'

Printed on 22 Feb 2024