

Datasheet: PPP025

Description:	RECOMBINANT PIG INTERLEUKIN-10
Name:	IL-10
Other names:	CYTOKINE STIMULATING INHIBITORY FACTOR
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Pig
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose
External Database Links	UniProt: Q29055 Related reagents

Entrez Gene:[397106](#) IL10 [Related reagents](#)

Product Information	Recombinant Pig interleukin-10 (IL-10) , also known as cytokine synthesis inhibitory factor (CSIF), is a secreted anti inflammatory cytokine produced primarily by monocytes. This product is a purified recombinant cytokine corresponding to aa 18-175 of porcine IL-10. IL-10 inhibits the synthesis of pro-inflammatory cytokines such as IFN γ , IL-2, IL-4, GM-CSF and TNF α .
Protein Molecular Weight	18.1 kDa
Purity	>95% by SDS PAGE analysis
ELISA	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2393 as a coating antibody and AHP2393B as a detection antibody.
Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10591 available at: 10591: https://www.bio-rad-antibodies.com/uploads/MSDS/10591.pdf
Regulatory	For research purposes only

North & South Tel: +1 800 265 7376**America** Fax: +1 919 878 3751Email: antibody_sales_us@bio-rad.com**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)

'M334929:181203'

Printed on 29 Aug 2021

© 2021 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)