

Datasheet: PPP032

**BATCH NUMBER 153956**

<b>Description:</b>	RECOMBINANT PIG INTERLEUKIN-3
<b>Name:</b>	IL-3
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Pig
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose
<b>Product Information</b>	<b>Recombinant Pig Interleukin-3 (IL-3)</b> is a recombinant protein corresponding to aa 24-144 of native pig IL-3.

IL-3 stimulates proliferation of granulocytes, monocytes, and dendritic cells in combination with other cytokines such as GM-CSF.

---

<b>Protein Molecular Weight</b>	13.7 kDa
---------------------------------	----------

---

<b>Purity</b>	>95% by SDS PAGE analysis
---------------	---------------------------

---

<b>ELISA</b>	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with <a href="#">AHP2819</a> as a coating antibody and <a href="#">AHP2819B</a> as a detection antibody.
--------------	--

---

<b>Storage</b>	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
----------------	--

---

<b>Guarantee</b>	12 months from date of despatch
------------------	---------------------------------

---

<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10591 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PPP032">https://www.bio-rad-antibodies.com/SDS/PPP032</a> 10591
--------------------------------------	---

---

<b>Regulatory</b>	For research purposes only
-------------------	----------------------------

---

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)  
'M373620:200923'

Printed on 19 Jan 2024

---

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)