

Datasheet: PRABP002

BATCH NUMBER 172906

Description:	RECOMBINANT RABBIT INTERFERON GAMMA
Name:	IFN GAMMA
Other names:	IFNG
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Source	<i>Pichia pastoris</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

**External Database
Links**

UniProt:

[P30123](#) [Related reagents](#)

Entrez Gene:

[100008602](#) IFNG [Related reagents](#)

Product Information **Recombinant Rabbit interferon gamma**, (IFN γ), known previously as macrophage activating factor, belongs to the type II interferon family and is produced by T lymphocytes and NK cells.

Native rabbit interferon gamma is a 167 amino acid cytokine with potent activating, anti-viral and anti-proliferative properties.

IFN γ is predominantly secreted by activated T lymphocytes in response to specific mitogens as a result of infection and has numerous immunological regulatory functions. Properties include macrophage and NK cell activation, upregulation of MHC class I and II molecules on various cells and induction of cytotoxicity.

**Protein Molecular
Weight**

16.9 kDa

Purity

>95% by SDS PAGE analysis

ELISA

This reagent may be used as a standard in a sandwich ELISA assay.

Storage

Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10591 available at:
<https://www.bio-rad-antibodies.com/SDS/PRABP002>

Regulatory

For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M426661:240213'

Printed on 29 Jan 2026