

Datasheet: PRABP006

BATCH NUMBER 166636

| | |
|----------------------|----------------------------------|
| Description: | RECOMBINANT RABBIT INTERLEUKIN-8 |
| Name: | IL-8 |
| Other names: | CXCL8 |
| Format: | Rec. Protein |
| Product Type: | Recombinant Protein |
| Quantity: | 5 µg |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|------------------|-----|----|----------------|--------------------|
| ELISA | ▪ | | | |
| Western Blotting | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

| | |
|---------------------------------|--|
| Target Species | Rabbit |
| Product Form | Purified recombinant protein - lyophilized |
| Reconstitution | Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution. |
| Preparation | Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography |
| Source | <i>Pichia pastoris</i> |
| Buffer Solution | Phosphate buffered saline |
| Preservative Stabilisers | 10% Trehalose |

**External Database
Links**

UniProt:

[P19874](#) [Related reagents](#)

Entrez Gene:

[100009129](#) IL8 [Related reagents](#)

Synonyms CXCL8

Product Information **Recombinant Rabbit interleukin-8 (IL-8)** is a recombinant protein corresponding to aa 23-101 of native rabbit IL-8.

Rabbit IL-8, also known as CXCL8, is a chemokine secreted by a range of cell types including macrophages and endothelial cells in response to infection and functions as a chemotactin on T lymphocytes, neutrophils and basophils.

Protein Molecular Weight 9.1 kDa

Purity >95% by SDS PAGE analysis

ELISA This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2992](#) as a coating antibody and [AHP2992B](#) as a detection antibody.

Storage Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

Guarantee 12 months from date of despatch

Health And Safety Information Material Safety Datasheet documentation #10591 available at: <https://www.bio-rad-antibodies.com/SDS/PRABP006>
10591

Regulatory For research purposes only

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M415639:230113'

Printed on 18 Jan 2024