

Datasheet: PRABP006

**BATCH NUMBER 169672**

<b>Description:</b>	RECOMBINANT RABBIT INTERLEUKIN-8
<b>Name:</b>	IL-8
<b>Other names:</b>	CXCL8
<b>Format:</b>	Rec. Protein
<b>Product Type:</b>	Recombinant Protein
<b>Quantity:</b>	5 µg

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Target Species</b>	Rabbit
<b>Product Form</b>	Purified recombinant protein - lyophilized
<b>Reconstitution</b>	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
<b>Preparation</b>	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
<b>Source</b>	<i>Pichia pastoris</i>
<b>Buffer Solution</b>	Phosphate buffered saline
<b>Preservative Stabilisers</b>	10% Trehalose

**External Database  
Links**

**UniProt:**

[P19874](#)    [Related reagents](#)

**Entrez Gene:**

[100009129](#)    IL8    [Related reagents](#)

---

**Synonyms**                    CXCL8

---

**Product Information**    **Recombinant Rabbit interleukin-8 (IL-8)** is a recombinant protein corresponding to aa 23-101 of native rabbit IL-8.

Rabbit IL-8, also known as CXCL8, is a chemokine secreted by a range of cell types including macrophages and endothelial cells in response to infection and functions as a chemotactin on T lymphocytes, neutrophils and basophils.

---

**Protein Molecular  
Weight**                    9.1 kDa

---

**Purity**                    >95% by SDS PAGE analysis

---

**ELISA**                    This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2992](#) as a coating antibody and [AHP2992B](#) as a detection antibody.

---

**Storage**                    Prior to reconstitution store at -20°C.  
After reconstitution store at -20°C.  
Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.

---

**Guarantee**                    12 months from date of despatch

---

**Health And Safety  
Information**                Material Safety Datasheet documentation #10591 available at:  
<https://www.bio-rad-antibodies.com/SDS/PRABP006>

---

**Regulatory**                    For research purposes only

---

**Product inquiries:** [www.bio-rad-antibodies.com/technical-support](http://www.bio-rad-antibodies.com/technical-support)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](http://bio-rad-antibodies.com/datasheets)

'M415639:230113'

**Printed on 29 Jan 2026**