

Datasheet: PRABP008

BATCH NUMBER 161844

Description:	RECOMBINANT RABBIT MCP-1
Name:	MCP-1
Other names:	CCL2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose
External Database Links	UniProt:

Entrez Gene:

[100009130](#) CCL2 [Related reagents](#)

Synonyms	SCYA2
Product Information	Recombinant Rabbit MCP-1 is an 11.4 kDa peptide corresponding to aa 24 – 125 of native rabbit monocyte chemotactic protein 1 (MCP-1). Also known as CCL2 , MCP-1 is a secreted chemokine involved in the recruitment of monocytes to sites of injury or inflammation. MCP-1 plays an important role in the inflammatory response of blood monocytes and tissue macrophages.
Protein Molecular Weight	11.4 kDa
Purity	>95% by SDS PAGE analysis
ELISA	This reagent may be used as a standard in a sandwich ELISA assay in conjunction with AHP2994 as a coating antibody and AHP2994B as a detection antibody.
Storage	Prior to reconstitution store at -20°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted.
Guarantee	12 months from date of despatch
Health And Safety Information	Material Safety Datasheet documentation #10591 available at: https://www.bio-rad-antibodies.com/SDS/PRABP008 10591
Regulatory	For research purposes only

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M334941:181203'

Printed on 18 Jan 2024