

Datasheet: PRABP008

BATCH NUMBER 170740

Description:	RECOMBINANT RABBIT MCP-1
Name:	MCP-1
Other names:	CCL2
Format:	Rec. Protein
Product Type:	Recombinant Protein
Quantity:	5 µg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
ELISA	▪			
Western Blotting	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Target Species	Rabbit
Product Form	Purified recombinant protein - lyophilized
Reconstitution	Reconstitute with 20 µl sterile PBS containing 0.1% carrier protein Care should be taken during reconstitution as the protein may appear as a film at the bottom of the vial. Bio-Rad recommend that the vial is gently mixed after reconstitution.
Preparation	Recombinant protein expressed in <i>Pichia pastoris</i> and purified by ion-exchange chromatography
Source	<i>Pichia pastoris</i>
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	10% Trehalose

**External Database
Links**

UniProt:

[P28292](#) [Related reagents](#)

Entrez Gene:

[100009130](#) CCL2 [Related reagents](#)

Synonyms

SCYA2

Product Information

Recombinant Rabbit MCP-1 is an 11.4 kDa peptide corresponding to aa 24 – 125 of native rabbit monocyte chemotactic protein 1 (MCP-1). Also known as CCL2 , MCP-1 is a secreted chemokine involved in the recruitment of monocytes to sites of injury or inflammation. MCP-1 plays an important role in the inflammatory response of blood monocytes and tissue macrophages.

**Protein Molecular
Weight**

11.4 kDa

Purity

>95% by SDS PAGE analysis

ELISA

This reagent may be used as a standard in a sandwich ELISA assay in conjunction with [AHP2994](#) as a coating antibody and [AHP2994B](#) as a detection antibody.

Storage

Prior to reconstitution store at -20°C.
After reconstitution store at -20°C.
Storage in frost-free freezers is not recommended. This product should be stored undiluted.

Guarantee

12 months from date of despatch

**Health And Safety
Information**

Material Safety Datasheet documentation #10591 available at:
<https://www.bio-rad-antibodies.com/SDS/PRABP008>

Regulatory

For research purposes only

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M415637:230113'

Printed on 29 Jan 2026
