

## Datasheet: PUR002

**BATCH NUMBER 21431**

<b>Description:</b>	PROTEUS PROTEIN A MIDI PURIFICATION SPIN COLUMN PACK
<b>Name:</b>	ANTIBODY PURIFICATION KIT PROTEIN A
<b>Format:</b>	MIDI - Spin Columns
<b>Product Type:</b>	Purification Kit
<b>Quantity:</b>	12 UNITS

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Preservative Stabilisers</b>	0.01% Thiomersal
---------------------------------	------------------

**Product Information** See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

<b>Reagents In The Kit</b>	12 x Protein A Spin Column Plugs 12 x Spin Columns 24 x 50 ml Centrifuge Tubes Plug Insertion Tool Instructions
----------------------------	---

**References** 1. Raeisi, H. *et al.* (2019) Applying the pthA effector protein of *Xanthomonas citri*. subsp. *citr.i* for production of specific antibodies and its application for detection of infected plants [Journal of Plant Pathology. 102 \(1\): 79-87.](#)

**Storage** Store at +4°C. DO NOT FREEZE.

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #qtec44-3 available at: qtec44-3: <a href="https://www.bio-rad-antibodies.com/uploads/MSDS/qtec44-3.pdf">https://www.bio-rad-antibodies.com/uploads/MSDS/qtec44-3.pdf</a> Not required
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN G MIDI PURIFICATION SPIN COLUMN PACK \(PUR011\)](#)

**North & South** Tel: +1 800 265 7376

**America** Fax: +1 919 878 3751

Email: [antibody\\_sales\\_us@bio-rad.com](mailto:antibody_sales_us@bio-rad.com)

**Worldwide**

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: [antibody\\_sales\\_uk@bio-rad.com](mailto:antibody_sales_uk@bio-rad.com)

**Europe**

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: [antibody\\_sales\\_de@bio-rad.com](mailto:antibody_sales_de@bio-rad.com)

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)

'M390557:210922'

**Printed on 09 Feb 2023**

---

© 2023 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)