

Datasheet: PUR002

BATCH NUMBER 230507

Description:	PROTEUS PROTEIN A MIDI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN A
Format:	MIDI - Spin Columns
Product Type:	Purification Kit
Quantity:	12 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers	0.01% Thiomersal
---------------------------------	------------------

Product Information	See Binding affinities for Protein A and G
----------------------------	--

Recommended Reading - [Proteus Protein A Antibody Purification Handbook](#)

Reagents In The Kit	12 x Protein A Spin Column Plugs 12 x Spin Columns 24 x 50 ml Centrifuge Tubes Plug Insertion Tool Instructions
----------------------------	---

References	1. Raeisi, H. <i>et al.</i> (2019) Applying the pthA effector protein of <i>Xanthomonas citri</i> . subsp. <i>citr.i</i> for production of specific antibodies and its application for detection of infected plants Journal of Plant Pathology. 102 (1): 79-87.
-------------------	---

Storage	Store at +4°C. DO NOT FREEZE.
----------------	-------------------------------

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information Material Safety Datasheet documentation #qtec44-3 available at:
<https://www.bio-rad-antibodies.com/SDS/PUR002>
qtec44-3
Not required

Regulatory For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN G MIDI PURIFICATION SPIN COLUMN PACK \(PUR011\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350455:190307'

Printed on 29 Feb 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)