

Datasheet: PUR002

BATCH NUMBER 25281

Description:	PROTEUS PROTEIN A MIDI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN A
Format:	MIDI - Spin Columns
Product Type:	Purification Kit
Quantity:	12 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

20% Ethanol

Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit

12 x Protein A Spin Column Plugs
 12 x Spin Columns
 24 x 50 ml Centrifuge Tubes
 Plug Insertion Tool
 Instructions

References

1. Raeisi, H. *et al.* (2019) Applying the pthA effector protein of *Xanthomonas citri*. subsp. *citri*.i for production of specific antibodies and its application for detection of infected plants [Journal of Plant Pathology. 102 \(1\): 79-87.](#)

Storage

This product is shipped at ambient temperature.

Store at +4°C. DO NOT FREEZE.

Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #20510 available at: https://www.bio-rad-antibodies.com/SDS/PUR002
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M441571:250523'

Printed on 08 May 2026

© 2026 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)