

Datasheet: PUR007

BATCH NUMBER 20609

Description:	PROTEUS PROTEIN A MINI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN A
Format:	MINI - Spin Columns
Product Type:	Purification Kit
Quantity:	48 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers	0.01% Thiomersal
---------------------------------	------------------

Product Information See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit	48 x Protein A Spin Column Plugs 48 x Spin Columns 96 x 2.2 ml Centrifuge Tubes Plug Insertion Tool Instructions
----------------------------	--

References

1. Bekesova, S. *et al.* (2012) N-glycans in liver-secreted and immunoglobulin-derived protein fractions. [J Proteomics. 75 \(7\): 2216-24.](#)
2. Goldman, R. (2012) Glycomic Analysis of Prostate Cancer Medicine and Medical Research Award Number: W81XWH-09-1-0202
3. Raeisi, H. *et al.* (2019) Applying the pthA effector protein of *Xanthomonas citri*. subsp.

citri. for production of specific antibodies and its application for detection of infected plants
[Journal of Plant Pathology. 102 \(1\): 79-87.](#)

Storage	Store at +4°C. DO NOT FREEZE.
----------------	-------------------------------

Guarantee	Guaranteed until date of expiry. Please see product label.
------------------	--

Health And Safety Information	Material Safety Datasheet documentation #qtec44-4 available at: https://www.bio-rad-antibodies.com/SDS/PUR007 qtec44-4 Not required
--------------------------------------	--

Regulatory	For research purposes only
-------------------	----------------------------

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN G MINI PURIFICATION SPIN COLUMN PACK \(PUR015\)](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M350458:190307'

Printed on 18 Oct 2024

© 2024 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)