

Datasheet: PUR007

BATCH NUMBER 24441

Description:	PROTEUS PROTEIN A MINI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN A
Format:	MINI - Spin Columns
Product Type:	Purification Kit
Quantity:	48 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

20% Ethanol

Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit

48 x Protein A Spin Column Plugs
 48 x Spin Columns
 96 x 2.2 ml Centrifuge Tubes
 Plug Insertion Tool
 Instructions

References

1. Bekesova, S. *et al.* (2012) N-glycans in liver-secreted and immunoglobulin-derived protein fractions. [J Proteomics. 75 \(7\): 2216-24.](#)
2. Raeisi, H. *et al.* (2019) Applying the pthA effector protein of *Xanthomonas citri*. subsp. *citri*. for production of specific antibodies and its application for detection of infected plants [Journal of Plant Pathology. 102 \(1\): 79-87.](#)

Storage	Store at +4°C. DO NOT FREEZE.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #20510 available at: https://www.bio-rad-antibodies.com/SDS/PUR007
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN G MINI PURIFICATION SPIN COLUMN PACK \(PUR015\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M433389:241018'

Printed on 08 May 2026