

## Datasheet: PUR011

**BATCH NUMBER 21221**

<b>Description:</b>	PROTEUS PROTEIN G MIDI PURIFICATION SPIN COLUMN PACK
<b>Name:</b>	ANTIBODY PURIFICATION KIT PROTEIN G
<b>Format:</b>	MIDI - Spin Columns
<b>Product Type:</b>	Purification Kit
<b>Quantity:</b>	12 UNITS

### Product Details

#### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

#### Preservative Stabilisers

20% Ethanol

#### Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

#### Reagents In The Kit

12 x Protein G Spin Column Plugs  
 12 x Spin Columns  
 24 x 50 ml Centrifuge Tubes  
 Plug Insertion Tool  
 Instructions

#### References

1. Montaudon, E. *et al.* (2015) Cardiac effects of long-term active immunization with the second extracellular loop of human  $\beta$ 1- and/or  $\beta$ 3-adrenoceptors in Lewis rats. [Pharmacol Res. 100: 210-9.](#)

#### Storage

Store at +4°C. DO NOT FREEZE

---

<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #qtec44-5 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PUR011qtec44-5">https://www.bio-rad-antibodies.com/SDS/PUR011qtec44-5</a> Not required
<b>Regulatory</b>	For research purposes only

---

## Related Products

### Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN A MIDI PURIFICATION SPIN COLUMN PACK \(PUR002\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	---	------------------	---	---------------	---

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
'M350460:190307'

**Printed on 18 Jan 2024**