

Datasheet: PUR014

BATCH NUMBER 210323

Description:	PROTEUS PROTEIN G MINI PURIFICATION STARTER KIT
Name:	ANTIBODY PURIFICATION KIT PROTEIN G
Format:	MINI - Spin Columns, Buffers, U/F Spinners
Product Type:	Purification Kit
Quantity:	2 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers	Buffers contain 0.09% Sodium Azide (NaN ₃) Resin contains 20% Ethanol
---------------------------------	--------------------------------------------------------------------------------------

Product Information See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein G Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit	2 x Protein G Spin Column Plugs 2 x Spin Columns 4 x 2.2 ml Centrifuge Tubes 2 x 10 kDa MWCO Ultrafiltration Spinners 1 x 60 ml of 1 x Buffer G 1 x 20 ml of 1 x Buffer B2 1 x 10 ml of 1 x Buffer C Plug Insertion Tool Instructions
----------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Storage	Store at +4°C. DO NOT FREEZE.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #10125 #10253 #10254 available at: https://www.bio-rad-antibodies.com/SDS/PUR014 Protein G Binding Buffer G (10125) Elution Buffer B2 (10253) Neutralisation Buffer C (10254)
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

- [MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)
- [RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)
- [PROTEUS PROTEIN A MINI PURIFICATION STARTER KIT \(PUR006\)](#)
- [MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets
'M350462:190307'

Printed on 18 Jan 2024