

Datasheet: PUR014

**BATCH NUMBER 220108**

<b>Description:</b>	PROTEUS PROTEIN G MINI PURIFICATION STARTER KIT
<b>Name:</b>	ANTIBODY PURIFICATION KIT PROTEIN G
<b>Format:</b>	MINI - Spin Columns, Buffers, U/F Spinners
<b>Product Type:</b>	Purification Kit
<b>Quantity:</b>	2 UNITS

## Product Details

### Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit [www.bio-rad-antibodies.com/protocols](http://www.bio-rad-antibodies.com/protocols).

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

<b>Preservative Stabilisers</b>	Buffers contain 0.09% Sodium Azide (NaN <sub>3</sub> ) Resin contains 20% Ethanol
---------------------------------	--------------------------------------------------------------------------------------

**Product Information** See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein G Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

<b>Reagents In The Kit</b>	2 x Protein G Spin Column Plugs 2 x Spin Columns 4 x 2.2 ml Centrifuge Tubes 2 x 10 kDa MWCO Ultrafiltration Spinners 1 x 60 ml of 1 x Buffer G 1 x 20 ml of 1 x Buffer B2 1 x 10 ml of 1 x Buffer C Plug Insertion Tool Instructions
----------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

<b>Storage</b>	Store at +4°C. DO NOT FREEZE.
<b>Guarantee</b>	Guaranteed until date of expiry. Please see product label.
<b>Health And Safety Information</b>	Material Safety Datasheet documentation #10125 #10253 #10254 available at: <a href="https://www.bio-rad-antibodies.com/SDS/PUR014">https://www.bio-rad-antibodies.com/SDS/PUR014</a> Protein G Binding Buffer G (10125) Elution Buffer B2 (10253) Neutralisation Buffer C (10254)
<b>Regulatory</b>	For research purposes only

## Related Products

### Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)  
[PROTEUS PROTEIN A MINI PURIFICATION STARTER KIT \(PUR006\)](#)  
[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

<b>North &amp; South America</b>	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: <a href="mailto:antibody_sales_us@bio-rad.com">antibody_sales_us@bio-rad.com</a>	<b>Worldwide</b>	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: <a href="mailto:antibody_sales_uk@bio-rad.com">antibody_sales_uk@bio-rad.com</a>	<b>Europe</b>	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: <a href="mailto:antibody_sales_de@bio-rad.com">antibody_sales_de@bio-rad.com</a>
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://bio-rad-antibodies.com/datasheets)  
 'M350462:190307'

Printed on 13 Aug 2023