

Datasheet: PUR015

BATCH NUMBER 24211

| | |
|----------------------|--|
| Description: | PROTEUS PROTEIN G MINI PURIFICATION SPIN COLUMN PACK |
| Name: | ANTIBODY PURIFICATION KIT PROTEIN G |
| Format: | MINI - Spin Columns |
| Product Type: | Purification Kit |
| Quantity: | 48 UNITS |

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

| | Yes | No | Not Determined | Suggested Dilution |
|----------------------|-----|----|----------------|--------------------|
| Protein Purification | ▪ | | | |

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

20% Ethanol

Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein G Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit

48 x Protein G Spin Column Plugs
 48 x Spin Columns
 96 x 2.2 ml Centrifuge Tubes
 Plug Insertion Tool
 Instructions

References

1. Bekesova S *et al.* (2012) N-glycans in liver-secreted and immunoglobulin-derived protein fractions. [J Proteomics. 75 \(7\): 2216-24.](#)
2. Goldman, R. (2012) Glycomic Analysis of Prostate Cancer. Medicine and Medical Research Award Number: W81XWH-09-1-0202

| | |
|--------------------------------------|--|
| Storage | Store at +4°C. DO NOT FREEZE. |
| Guarantee | Guaranteed until date of expiry. Please see product label. |
| Health And Safety Information | Material Safety Datasheet documentation #20510 available at: https://www.bio-rad-antibodies.com/SDS/PUR015 20510 Not required |
| Regulatory | For research purposes only |

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN A MINI PURIFICATION SPIN COLUMN PACK \(PUR007\)](#)

| | | | | | |
|----------------------------------|---|------------------|---|---------------|---|
| North & South America | Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com | Worldwide | Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com | Europe | Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com |
|----------------------------------|---|------------------|---|---------------|---|

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M417829:230417'

Printed on 18 Oct 2024