

Datasheet: PUR015

BATCH NUMBER 25491

Description:	PROTEUS PROTEIN G MINI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN G
Format:	MINI - Spin Columns
Product Type:	Purification Kit
Quantity:	48 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

20% Ethanol

Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein G Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit

48 x Protein G Spin Column Plugs
 48 x Spin Columns
 96 x 2.2 ml Centrifuge Tubes
 Plug Insertion Tool
 Instructions

References

1. Bekesova S *et al.* (2012) N-glycans in liver-secreted and immunoglobulin-derived protein fractions. [J Proteomics. 75 \(7\): 2216-24.](#)
2. Goldman, R. (2012) Glycomic Analysis of Prostate Cancer. Medicine and Medical Research Award Number: W81XWH-09-1-0202

Storage	This product is shipped at ambient temperature. Store at +4°C. DO NOT FREEZE.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #20510 available at: https://www.bio-rad-antibodies.com/SDS/PUR015
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN A MINI PURIFICATION SPIN COLUMN PACK \(PUR007\)](#)

Product inquiries: www.bio-rad-antibodies.com/technical-support

To find a batch/lot specific datasheet for this product, please use our online search tool at: bio-rad-antibodies.com/datasheets

'M441579:250523'

Printed on 08 May 2026