

Datasheet: PUR015

BATCH NUMBER 310712

Description:	PROTEUS PROTEIN G MINI PURIFICATION SPIN COLUMN PACK
Name:	ANTIBODY PURIFICATION KIT PROTEIN G
Format:	MINI - Spin Columns
Product Type:	Purification Kit
Quantity:	48 UNITS

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Protein Purification	▪			

Where this product has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the product for use in their own system using appropriate negative/positive controls.

Preservative Stabilisers

20% Ethanol

Product Information

See [Binding affinities for Protein A and G](#)

Recommended Reading - [Proteus Protein A Purification Handbook](#)

For more information visit [Proteus: Protein Purification Spin Kits](#)

Reagents In The Kit

48 x Protein G Spin Column Plugs
48 x Spin Columns
96 x 2.2 ml Centrifuge Tubes
Plug Insertion Tool
Instructions

References

1. Bekesova S *et al.* (2012) N-glycans in liver-secreted and immunoglobulin-derived protein fractions. [J Proteomics. 75 \(7\): 2216-24.](#)
2. Goldman, R. (2012) Glycomic Analysis of Prostate Cancer. Medicine and Medical Research Award Number: W81XWH-09-1-0202

Storage	Store at +4°C. DO NOT FREEZE.
Guarantee	Guaranteed until date of expiry. Please see product label.
Health And Safety Information	Material Safety Datasheet documentation #qtec44-6 available at: https://www.bio-rad-antibodies.com/SDS/PUR015qtec44-6 Not required
Regulatory	For research purposes only

Related Products

Recommended Useful Reagents

[MOUSE MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(MMT1\)](#)

[RAT MONOCLONAL ANTIBODY ISOTYPING TEST KIT \(RMT1\)](#)

[MYCOPLASMA REMOVAL AGENT \(BUF035\)](#)

[PROTEUS PROTEIN A MINI PURIFICATION SPIN COLUMN PACK \(PUR007\)](#)

North & South America	Tel: +1 800 265 7376 Fax: +1 919 878 3751 Email: antibody_sales_us@bio-rad.com	Worldwide	Tel: +44 (0)1865 852 700 Fax: +44 (0)1865 852 739 Email: antibody_sales_uk@bio-rad.com	Europe	Tel: +49 (0) 89 8090 95 21 Fax: +49 (0) 89 8090 95 50 Email: antibody_sales_de@bio-rad.com
----------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------	------------------	-------------------------------------------------------------------------------------------------------------------------------------------------	---------------	-----------------------------------------------------------------------------------------------------------------------------------------------------

To find a batch/lot specific datasheet for this product, please use our online search tool at: [bio-rad-antibodies.com/datasheets](https://www.bio-rad-antibodies.com/datasheets)
'M350463:190307'

Printed on 18 Jan 2024