

Datasheet: STAR116

Description:	GOAT ANTI RAT IgM
Specificity:	IgM
Format:	Purified
Product Type:	Polyclonal Antibody
Isotype:	Polyclonal IgG
Quantity:	0.5 mg

Product Details

Applications

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, peer-reviewed publications or personal communications from the originators. Please refer to references indicated for further information. For general protocol recommendations, please visit www.bio-rad-antibodies.com/protocols.

	Yes	No	Not Determined	Suggested Dilution
Flow Cytometry	▪			
Immunohistology - Frozen			▪	
Immunohistology - Paraffin			▪	
ELISA	▪			1ug/ml - 10ug/ml
Immunoprecipitation			▪	
Western Blotting	▪			

Where this antibody has not been tested for use in a particular technique this does not necessarily exclude its use in such procedures. Suggested working dilutions are given as a guide only. It is recommended that the user titrates the antibody for use in their own system using appropriate negative/positive controls.

Target Species	Rat
Product Form	Purified IgG - liquid
Antiserum Preparation	Antisera to rat IgM were raised by repeated immunisations of goats with highly purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Buffer Solution	Phosphate buffered saline
Preservative Stabilisers	0.09% Sodium Azide
Approx. Protein Concentrations	IgG concentration 1.0 mg/ml
RRID	AB_567061
Specificity	Goat anti Rat IgM antibody recognizes rat IgM as assessed by immunoelectrophoresis and ELISA. Goat anti Rat IgM antibody may cross react with IgM from other species.

Flow Cytometry Use 10ul of the suggested working dilution to stain 1 x 10⁶ cells in 100ul.

References

1. Terrell, S.P. *et al.* (2012) Glomerulonephropathy in aged captive Key Largo woodrats (*Neotoma floridana smalli*). [Vet Pathol. 49 \(4\): 710-6.](#)
2. Kesik-brodacka, M. *et al.* (2017) Immune response of rats vaccinated orally with various plant-expressed recombinant cysteine proteinase constructs when challenged with *Fasciola hepatica* metacercariae. [PLoS Negl Trop Dis. 11 \(3\): e0005451.](#)
3. Bouma, H.R. *et al.* (2013) Hibernation is associated with depression of T-cell independent humoral immune responses in the 13-lined ground squirrel. [Dev Comp Immunol. 39 \(3\): 154-60.](#)

Storage Store at +4°C. DO NOT FREEZE.
This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Guarantee 12 months from date of despatch.

Health And Safety Information Material Safety Datasheet documentation #10040 available at:
10040: <https://www.bio-rad-antibodies.com/uploads/MSDS/10040.pdf>

Regulatory For research purposes only

Related Products

Recommended Secondary Antibodies

Rabbit Anti Goat IgG (Fc) (STAR122...) [FITC](#), [HRP](#)

North & South Tel: +1 800 265 7376

America Fax: +1 919 878 3751

Email: antibody_sales_us@bio-rad.com

Worldwide

Tel: +44 (0)1865 852 700

Fax: +44 (0)1865 852 739

Email: antibody_sales_uk@bio-rad.com

Europe

Tel: +49 (0) 89 8090 95 21

Fax: +49 (0) 89 8090 95 50

Email: antibody_sales_de@bio-rad.com

'M339611:181219'

Printed on 15 Oct 2019

© 2019 Bio-Rad Laboratories Inc | [Legal](#) | [Imprint](#)